# PrinterNOx

Operating Manual – English

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## Contents

I.	Important safety information – English .....	3
	Introduction .....	4
	Medical Use .....	5
	NO <sub>2</sub> Production .....	6
	Environmental Safety .....	6
	Operation .....	7
	1. Measurement .....	10
	2. Calibration .....	17
	3. Set Alarms .....	23
	4. Report to printer .....	24
	5. Status .....	24
	6. Setup .....	25
	7. Calculate Flow .....	26
	Battery operation .....	27
	Paper Loading .....	28
	Battery Maintenance .....	29
	Cleaning Procedure .....	31
	Servicing .....	31
	Symbols .....	32
	Electrical Environment .....	32
	Specifications .....	32
II.	Wichtige Sicherheitshinweise – Deutsch .....	36
	Einleitung .....	37
	Medizinische Anwendung .....	38
	NO <sub>2</sub> Produktion .....	39
	Umweltschutz .....	40
	Betrieb .....	40
	1. Messung .....	44
	2. Kalibration .....	52
	3. Alarmeinstellungen .....	58
	4. Bericht an den Drucker .....	59

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# PrinterNOx

## Operating manual

### **I. Important safety information - English**

Nitric oxide is a highly toxic gas, exposure to levels as low as 200 parts per million (ppm) may be fatal even after very short exposures, and prolonged exposure to much lower levels may be harmful. All users of nitric oxide must familiarise themselves with the appropriate safety regulations. When calibrating the PrinterNOx monitor it is important to ensure that calibration gas is adequately vented and safely exhausted. Nitric oxide reacts with oxygen in air to form nitrogen dioxide, an extremely toxic gas with a maximum allowable industrial ambient level of 5ppm.

The PrinterNOx nitric oxide/nitrogen dioxide monitor has been designed to the highest standards and incorporates many inbuilt self-checks. As with any measurement device, however, malfunction is possible. When used in the administration of nitric oxide the PrinterNOx must not, therefore, be incorporated as a primary part of the NO delivery system, but rather as a safety check on the delivery system. The delivery system should be constructed in such a way as to minimise possible errors. CareFusion Ltd does not accept any liability for injury resulting from malfunction and failure to follow these guidelines. Failure to follow instructions in the remainder of this manual may also result in dangerously inaccurate measurements being obtained.

## Introduction

The PrinterNOx is a combined Nitric Oxide / Nitrogen Dioxide monitor based on electrochemical fuel cells and may be operated directly from the mains power or from the internal rechargeable NiCad battery pack. Fuel cells work through the reaction of the target gas with an electrolyte at one electrode and oxygen (from ambient air) at the other. This reaction generates an electrical current proportional to gas concentration. In comparison with other technologies (such as chemiluminescent detectors) electrochemical sensors are compact, inexpensive and much easier to use with zero warm up time and simple calibration. The PrinterNOx is designed to monitor delivered NO used as a therapy for adult respiratory distress syndrome and pulmonary hypertension on ventilated patients. The more toxic and non therapeutic by-product, NO<sub>2</sub>, is also monitored. The NO measurement has a resolution of 0.05 parts per million (ppm) with a full scale of 100ppm and the NO<sub>2</sub> measurement has a resolution of 0.05ppm with a full scale of 20ppm. The gas to be monitored is sampled using a constant flow (250ml/min), low dead space, sidestream system. The system includes a water trap, with integral hydrophobic filter, and humidity conditioning circuit to protect the fuel cells against extremes of humidity and particulate contamination. The measurements are displayed on a large graphic liquid crystal display and can be printed as a chart, on the integral thermal printer, together with calibration due date and alarm settings. The last 24 hours of NO and NO<sub>2</sub> measurements taken are recorded and stored in the internal battery backed memory and can be printed or uploaded on the RS232 interface for computer storage. The PrinterNOx includes a number of safety features and alarms. These include upper and lower alarms for NO and an upper alarm for NO<sub>2</sub>. An external alarm cable is available, and when connected to a NO delivery system, is used to trigger an alarm or initiate a safe operating mode if high concentrations of gas are detected. This concentration is 100ppm for NO, 9ppm NO<sub>2</sub> or alarm level whichever

is the greater. The flow of gas in the sampling line is continuously monitored and an alarm sounded if the line is occluded or the water trap is full. A low battery condition is also alarmed.

The PrinterNOx has been designed with ease of use and reliability as of foremost importance. To ensure accurate results, however, it is essential to familiarise yourself with the contents of this short manual.

## Medical Use

Nitric oxide (NO) is evolving as a novel therapy used to reduce pulmonary arterial pressure and improve gas exchange in paediatric and adult intensive care settings. The NO is usually added to the ventilator circuit from a cylinder of NO in nitrogen with a concentration up to 1000ppm and diluted in the circuit to between 1 and 40ppm.

Therapeutic concentrations must be well controlled as levels as low as 100ppm may cause lung injury.

For safety reasons sources with greater than 1000ppm NO should not be used. Lower concentrations carry less risk either to the patient in the case of accidental overdose, or to the hospital staff in the case of leakage. However, if moderately high doses of NO are to be administered then the decrease in the maximum FiO<sub>2</sub> must be considered and is given by the relationship:

$$\text{FiO}_2 (\text{max.}) = (1 - \text{Required NO Conc.} / \text{Source NO Conc.}) \times 100\%$$

Therefore, if 40ppm is required from a source of 500ppm NO in nitrogen then the maximum oxygen available would be:

$$\text{FiO}_2 (\text{max.}) = (1 - 40 / 500) \times 100\% = 92\%$$

Methods for therapeutic administration of nitric oxide depend on the exact type of ventilator and ventilator circuit to be used and have not yet been standardised. It is therefore the responsibility of the user to pay due regard to current medical practice when using the unit to monitor concentrations of inhaled nitric oxide and nitrogen dioxide.

## **NO<sub>2</sub> Production**

Nitric oxide in the presence of oxygen undergoes the following chemical reaction to form nitrogen dioxide, (a more toxic gas than nitric oxide):

The rate of conversion is highly dependent on the NO concentration and this has some important implications for the NO delivery system.

Firstly, sources of NO with excessively high concentrations should be avoided.

Secondly, the administered concentration should be kept as low as is clinically acceptable.

Finally, the administration system should mix the incoming NO quickly and effectively to rapidly reduce its concentration to the therapeutic level where NO<sub>2</sub> production will proceed at a much slower rate.

## **Environmental Safety**

To eliminate environmental pollution with NO and NO<sub>2</sub> the PrinterNOx is supplied with two NOxSORB chemical scavengers for removing these gases from the exhaust. The NOxSORB works by converting NO to NO<sub>2</sub>, which is then chemically absorbed.


Connect a NOxSORB to the exhaust of the PrinterNOx, as shown overleaf, orientated such that the arrow is pointing away from the instrument. The NOxSORB should be replaced after 24 hours or when the mauve spheres turn dark brown. When re-ordering please quote 36-SORB100 for a box of 10 or 36-SORB200 for a box of 20.

Environmental pollution with NO and NO<sub>2</sub> may also arise when leaks are present in the connections to the source of NO.

## Operation


Before operating the PrinterNOx, ensure that the unit has been fully charged using the AC adapter supplied. (See Charging Procedure.) The PrinterNOx may also be operated with the internal rechargeable batteries.

Connect the water trap with integral hydrophobic PTFE filter and sampling line to the PrinterNOx as shown below:


Switch on the PrinterNOx. **Please note that the water trap must be used at all times during measurement to avoid contamination of the fuel cells.**

The on/off switch is located at the rear of the unit:

An introductory screen giving the software version will be displayed briefly:

During this time the following information will be printed:


```
12/04/95 16:14  
NO calibration due : 01/01/96  
NO2 calibration due : 07/01/96  
NO upper alarm : 15  
NO lower alarm : 5  
NO2 upper alarm : 4
```

This information can be checked and the date, time, calibration or alarm settings may be adjusted prior to performing a measurement.


(See sections on **Setup, Calibration, and Set Alarms**)

The calibration due date is set to occur thirty days after the last calibration. However, calibration may be performed at any time before this date to improve measurement accuracy. Measurements can be made after the calibration due date has passed but consideration must be given to the gas sensors which exhibit a slow decline in sensitivity over time. This loss of sensitivity is typically less than 2% per month but could be more if they are exposed to adverse conditions such as exposure to gas concentrations above the measurement range or to solvents such as alcohol.

If the calibration due date has passed then the following reminder will be displayed both after the introductory screen and when commencing measurement:

Once ENTER is pressed the first menu is displayed and allows selection of several options:

Whilst this display is active the sampling pump is turned off, measurements are not taken, and the gas level alarms are inactive.

If no options are selected then the unit defaults to option 1 "Measurement" after 45 seconds. This is done so that the unit defaults to monitoring gas levels with alarms activated in the event that the instrument was inadvertently left on the first screen or was switched off and on.

During all other operations, a short warning alarm will occur, after 45 seconds of inactivity, to remind the user that measurement mode is not operational.

## **1. Measurement**

To obtain maximum accuracy from the measurement careful consideration must be given to the NO administration system in use.


Gas measurement with ventilators producing cyclic flows may give an overestimate or under estimate of the average NO concentration.

This is due to localised variation of NO concentration within the delivery system. Dependent on the measurement point, therefore, the sampled average may not correspond to the delivered average. It is therefore important to ensure that the NO delivery system is designed to avoid large variations in the level of delivered NO.

The PrinterNOx utilises a sidestream gas system, which is servo controlled to sample at a constant rate of 250 ml/min. By supplying the gas sensors with this constant flow, a high degree of accuracy and sensitivity can be achieved.


To carry out NO / NO<sub>2</sub> measurement, press "1".

The sampling pump will start and the following screen will be displayed:

The system should now be tested for leaks by occluding the sampling line. The pump will be heard to rise in pitch as the servo system tries to overcome the restriction.

As the pump cannot maintain the required flow rate then the measurement will cease to be displayed, the pump will continue to run, and the following message will be displayed:

When the occlusion is removed, the display will automatically revert to the measurement screen. Repeat this procedure by occluding the exhaust port.

The sampling line should now be connected to the administered gas, normally as close to the patient as possible. The sampling port must be positioned in such a way that water cannot enter the tubing, causing rapid filling of the water trap.

If the gas line appears to be continuously occluded, remove the water trap and empty any water into a suitable container. Ensure that there is no occlusion on the exhaust port.

Under exceptional circumstances, water may pass through a damaged water trap and obstruct the inlet to the PrinterNOx. In this case turn the unit off and draw any liquid out of the inlet with a syringe connected directly to the inlet. Replace the water trap with a new one, turn the unit on, and return to the measurement screen.

On the left-hand side of the NO legend are the upper and lower alarm levels for NO. On the left-hand side of the NO<sub>2</sub> legend is the upper alarm level for NO<sub>2</sub>.


On the bottom line, the rectangle will blink on or off every time a measurement is updated (every 2.6 seconds) to indicate correct functioning of the instrument.

The gas concentrations are given in parts per million (ppm).

If the printer option is turned on (see **Setup** section), then the following printout will be obtained during measurement:

```
12/04/95 16:14
NO calibration due : 01/01/96
NO2 calibration due : 07/01/96
NO upper alarm : 15
NO lower alarm : 5
NO2 upper alarm : 4
```

```
Patient :.....
Hospital :.....
Technician :.....
```

Spaces for entering the patient name or ID, hospital name, and technician name are provided.

This is followed by the axes for the NO and NO<sub>2</sub> graphs.

These graphs are auto scaled to the current alarm settings.

A section of the graph will be printed once every 32 samples and includes a time stamp every 15 minutes on the left-hand side. Each dot on the graph represents the average of 4 measurements.


If the measured gas goes outside the range of the alarm settings then an audible alarm and a red flashing visual alarm will be generated. In addition, the relevant indicator on the liquid crystal display will flash. The graph is then replaced by a printout of each individual measurement together with a time stamp and an alarm indicator as shown below:

```

12/04/95 16:14
NO calibration due : 01/01/96
NO2 calibration due : 07/01/96
NO upper alarm : 15
NO lower alarm : 5
NO2 upper alarm : 4
  
```

```

Patient :.....
Hospital :.....
Technician :.....
  
```


If the NO measurement exceeds 100ppm or NO<sub>2</sub> exceeds 9ppm then the external alarm, accessible at the rear of the instrument, is

activated. If the NO<sub>2</sub> alarm level is set above 9ppm, then the external alarm comes into operation at the set level for NO<sub>2</sub>. The external alarm consists of a pair of normally closed contacts, which are latched open when the above condition occurs. This can be used to sound an external alarm or switch a NO delivery system into safe operating mode in systems equipped with this facility.

The display will change to:


Measurements will continue but it should be noted that the measurement range would be limited to just over 100ppm NO and 20ppm NO<sub>2</sub> and that the actual concentrations could be significantly higher.

Press enter to clear the external alarm when the reading drops below the external alarm levels.

When this condition occurs, the cause should be immediately investigated. However, caution should be observed when shutting down the supply of NO, as the sudden cessation of NO can be more detrimental to the patient than a temporary overdose of NO.


Ideally, the external alarm should be used to switch the delivery system into a 'safe' operating mode or switch in a backup supply.

Should the water trap fill during the course of the measurement, it will overflow into the integral hydrophobic PTFE filter, cause an obstruction, and the following message will be displayed:

Unscrew the two retaining screws, and remove the water trap from the PrinterNOx.

Empty the accumulated water into a suitable container by connecting a syringe (20ml minimum size) with the adapter supplied and pushing air through the assembly as shown below.

Reconnect water trap to the PrinterNOx and secure with the two retaining screws. If the 'Gas Line Occluded' message reappears, disassemble and replace the water trap with a new one.


Under exceptional circumstances, water may pass through a damaged water trap and obstruct the inlet to the PrinterNOx. In this case turn the unit off and draw any liquid out of the inlet with a syringe connected directly to the inlet. Replace the water trap with a new one, turn the unit on, and return to the measurement screen.

The alarm may be muted for a 45-second period by pressing the CANCEL key.

## 2. Calibration


If the instrument is in use immediately before calibration is attempted, then the system should be purged with clean air. This may be done by disconnecting the water trap from the instrument and running the instrument in measurement mode for 5 minutes.

To carry out sensor calibration press 2 and the following screen will be displayed:

The PrinterNOx requires clean air free from any NO or NO<sub>2</sub> so that the sensors can be zeroed.


With the water trap disconnected, press enter and the display will change to:


The sampling pump will start and the PrinterNOx will monitor the output of both sensors until a stable reading is obtained.

During this time, a line of asterisks appears on the screen to indicate correct operation.

When zeroing is complete the readings for the NO and NO<sub>2</sub> sensors are written into the battery backed memory and the display is changed to:


Press "1" to calibrate NO and the display will change to:

Enter the required gas concentration using the DELETE key to correct any mistakes and press ENTER.


The display will change to:

CareFusion can supply calibration gas (25ppm NO in nitrogen and 10ppm NO<sub>2</sub> in air) complete with a control valve for simple and economic user calibration.


The gas is supplied in convenient, disposable, aluminium canisters containing 58 litres of gas (Cat No 36-GAS100 for 25ppm NO in nitrogen and 36-GAS200 for 10ppm NO<sub>2</sub> in air). A flow regulator pre-set to deliver 500ml/min is available separately, Cat No 36-GAS500. With careful use, approximately 45 calibrations can be obtained from each canister.

Screw the control valve on to the canister and connect to the PrinterNOx with the water trap connected as shown below:


Slowly turn the control knob fully anti-clockwise. This will then supply a flow of approximately 500ml/min. As the PrinterNOx samples at a rate of 250ml/min then 250ml/min will be vented through the one way valve.

Press ENTER on the keypad and the display will change to:

The sampling pump will start and the PrinterNOx will monitor the output of the NO sensor until a stable reading is obtained. When stability has been achieved the calibration factor, time and date are written into the battery backed memory and the display returns to the previous screen. Remove the control valve by turning anti-clockwise to stop the gas flow. The NO<sub>2</sub> sensor can now be calibrated using the same procedure if required.

Within the lifetime of the sensors, zero drift of more than 10ppm (1ppm for NO<sub>2</sub>) or a calibration drift of more than 40% implies abnormal ageing of the electrochemical cell. If these errors are detected during the zeroing procedure or the calibration adjustment then the following will be displayed:

The most likely reason for this message is an operational error in the calibration procedure rather than a faulty sensor.

If this occurs during zeroing repeat the procedure and check that there is no nitric oxide, or nitrogen dioxide contamination, and that the water trap has been removed.

If the calibration error message appears during sensitivity adjustment, repeat the calibration procedure and check:

1. Is calibration gas within its warranty period?
2. Has correct gas concentration been entered?
3. Have the correct connections been made with no leaks?
4. Is there sufficient calibration gas flow?
5. Is there any reaction between nitric oxide and the connectors used?
6. Are there any contaminants between nitric oxide supply and sensor?


Nitric oxide is a very reactive gas and it is difficult to maintain a stable concentration within a cylinder. It is therefore advisable to have some independent means of checking the calibration gas used.

After calibration is complete, selecting option "3 Exit" will return the display to the main menu screen.


### 3. Set Alarms

Selecting option three will allow the NO and NO<sub>2</sub> alarm levels to be set individually from the following menu:

The current alarm levels are displayed for reference and, for example, the upper alarm level for NO may be altered by selecting option 1 when the following screen will be displayed:

Type in the required alarm level and press ENTER.

The lower NO and NO<sub>2</sub> alarm levels are set with the same procedure with option 2 and 3 respectively.

The valid range of entry is 0 to 99 for the NO upper and lower alarms and the upper alarm must be set higher than the lower alarm.

For NO<sub>2</sub>, the valid range of entry is zero to 19

When the required alarms have been set, then selecting option 4 will return to the main menu.

#### **4. Report to printer (RS232)**

The last 24 hours of measurements are stored in battery backed memory and may be either printed or directed to the RS232 port at the rear of the instrument. The choice of recovery method depends upon whether the printer or RS232 port is currently active (see **Setup**).

If the RS232 port is active then the data, in ASCII format, may be uploaded to a PC running standard communication software. A serial cable, Cat No 36-CAB1108, is available from CareFusion.


Note: Please ensure that all PC's and peripheral equipment comply to EN60950 'Safety of Information Technology Equipment' and keep them out of reach of the patient when connected to the PrinterNOx.

#### **5. Status**

This option allows the current configuration and sensor information to be displayed.

## 6. Setup

Selecting this menu item allows various options to be set and displays the screen:

The first option clears the stored 24 hours data from the memory. This is done when the currently stored data is of no interest or has previously been accessed so that future access of required data will be quicker.

The second option is used to change the date and time stored in the real time clock.


“Toggle output” is used to select whether data is sent to the printer or the RS232 port at the rear of the instrument. An asterisk appears against the selected output device. The initial few lines of calibration and alarm information is always printed regardless whether the printer is selected.

The date format can be changed using option 4.

Option 5 returns the display to the main menu screen.


## 7. Calculate Flow

This facility is used to calculate the required NO flow to achieve the required therapeutic concentration for ventilators which employ a constant bias flow arrangement. When selected the following screen is displayed:

Enter ventilator  
bias flow:  
  
— L/min

The valid range for the bias flow is 1-99 and once entered the display changes to:

Enter NO gas  
concentration in  
bottle:  
  
— ppm

Gas concentrations of between 50 and 1000ppm may be entered.

Following this, the required NO concentration in the range 1 to 100 should be entered:

```
Enter required NO
gas concentration
to patient:

  _ ppm
```

When all the data is entered the required NO flow and the maximum FiO<sub>2</sub> is calculated and displayed on the following screen:


```
Bias flow      70 L/m
Bottle conc   650 ppm
Patient conc   20 ppm

Req NO flow    2.22 L/m
Maximum FiO2  96.9 %

Press key to continue
```

## Battery operation

During patient transportation or mains power failure, the PrinterNOx may be powered from the internal rechargeable batteries for approximately 2 hours.

When the batteries are low the alarm will sound continuously and the following screen will be displayed:

## **Paper Loading**

The PrinterNOx requires thermal paper (Part Number 36-PSA1600 - Pack of 5 rolls) for the printer.

A warning line appears when 30cm of the paper roll remains. A new roll should be loaded when this mark becomes visible.

To load paper:

- a. Switch on the PrinterNOx.
- b. Cut the leading edge of the paper so that it is parallel to the roll shaft.
- c. Remove the paper cover.
- d. Align the paper so that it is straight when it is inserted into the paper inlet.
- e. Press the "PAPER" key to feed the paper until the leading edge of the paper passes the paper cutter.
- f. Replace paper cover.

Make sure that the specified paper (Part Number 36-PSA1600) is used as excessive wear on the thermal print head may occur if paper from other manufacturers is used. The following precautions should be observed when storing or handling the thermal paper rolls:

- Do not allow the printer to run without paper loaded.
- Do not store paper rolls in places exposed to direct sunlight.
- Do not use pastes containing organic solvents such as alcohol, esters or ketones.
- Do not allow the paper to be exposed to vinyl chloride.
- Do not use cellophane tape on the paper.
- Do not allow printer paper to be exposed to diazo copy paper that has been processed.

## Battery Maintenance

The PrinterNOx is powered from a 7.2volt 600mA-hr battery pack located under the paper cover or from the AC adapter supplied.

The battery pack is fully discharged when shipped from the factory.

**Note:** Ensure that the dc plug is inserted before the mains supply is connected to the unit and that the main is switched off before removing the plug.

Please re-charge the pack with the AC adapter supplied before use.

The battery pack will exhibit a slow decrease in capacity at a rate dependant upon usage. If the PrinterNOx is normally powered from the AC adapter then the battery pack will be on continuous trickle charge. Under these conditions, the capacity will be greater than 60% of nominal capacity for a period of at least 5 years. If the PrinterNOx is normally powered from the battery pack then the capacity will be greater than 60% of nominal capacity for at least 1000 charge/discharge cycles. To obtain maximum life, the batteries

should be allowed to discharge until the battery low warning is activated and then fully charged as in the following procedure.

- a. Plug the adapter into the AC outlet and connect the cable to the printer.
- b. The Charging light at the front of the keypad will be illuminated while the batteries are being charged.
- c. The battery pack will take 16 hours to become fully charged.
- d. **Warning!** Use only the AC adapter supplied. Use of any other type may cause permanent damage to the PrinterNOx and may cause a fire or electrical hazard.

**Notes:**

- Do not leave the PrinterNOx on with discharged batteries for long periods since this may result in deterioration of the batteries.
- The room temperature should be between 5 to 40 degrees C during charging.
- The printer can be used immediately after the AC adapter is connected, however, if the battery pack is not fully charged, printing may stop temporarily until sufficient charge is available.
- Do not short-circuit the battery pack.
- Batteries at the end of their service life should be returned to CareFusion or an authorised representative for correct disposal. These must only be replaced with the approved part (CareFusion part No. 36-BAT1033)
- The PrinterNOx is CLASS I equipment (mains earth connected).
- **Warning!** Do not puncture or incinerate the rechargeable batteries.

The PrinterNOx also contains two primary lithium cells. One of these cells is associated with the gas sensors conditioning circuit and is replaced when the gas sensors expire in accordance with the service manual. If the PrinterNOx is to remain unused for a period greater than the service life of the gas sensors (18 months) then this battery and the gas sensors must be removed in accordance with the service manual. The other primary lithium cell is used for memory backup purposes with a design life of at least 10 years and must not be removed.

## **Cleaning Procedure**

The water trap with integral filter and sampling line are for single patient use only and must be discarded after use. The outside of the PrinterNOx may be cleaned using a light rubbing action with a damp cloth. The PrinterNOx must be removed from the mains supply during cleaning and should be completely dry before re-connection.

## **Looking after your PrinterNOx**

Please observe the following precautions:

- Do not keep the unit in a damp place or expose it to extremes of temperature.
- Do not expose any part of the gas sampling system to any solvents e.g. alcohol as these will damage the fuel cells.

## **Servicing**

A full service manual including circuit diagrams and parts list is available upon request. Servicing should only be carried out by suitably trained personnel.

## Symbols

Type B device

In accordance with Directive 93/42/EEC

Disposal in compliance with WEEE

## Electrical Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

In the event of interference it is advisable to remove the offending equipment or increase the distance. Failure to do so may compromise the safe operation of the PrinterNOx.

## Specifications

Gases detected		Nitric Oxide (NO)
		Nitrogen Dioxide (NO <sub>2</sub> )
Sensor type		Electrochemical fuel cells
Measurement ranges	NO:	100 Parts per million (ppm)
	NO <sub>2</sub> :	20ppm
Resolution	NO:	0.05ppm
	NO <sub>2</sub> :	0.05ppm

Accuracy	NO: +/-5% or 0.05ppm whichever is the greater. NO <sub>2</sub> : +/-5% or 0.05ppm whichever is the greater
Response time	NO: Less than 10 seconds (90% FSD) NO <sub>2</sub> : Less than 30 seconds (90% FSD)
Alarm range	NO: lower and upper 0 to 99ppm NO <sub>2</sub> : upper 0 to 19ppm
Safety cut-out operation	NO: 100ppm NO <sub>2</sub> : 9ppm
Alarmed conditions	Water trap full Sample line occluded. Battery low
Sensor life	>1 year
Sidestream flow rate	250 ml/minute
Input pressure range	0 to 100 cmH <sub>2</sub> O
Data storage	24 hours
Printer type	320 dot / line, low power thermal printer
Paper type	11cm thermal paper
Power supply (Class I Equipment)	Primary 100–250volt 50 - 60Hz Secondary 12volt 1.5A DC
Battery pack	Rechargeable NiCad 7.2V 600mA-hours
Dimensions	315 x 120 x 90 mm
Weight	1.2 Kg
Operating temperature:	0 to +40 °C
Operating humidity:	30% to 90% RH
Internal battery life	10 years
Operating pressure:	Atmospheric +/-10%
Storage temperature	-20° to +70 °C
Storage humidity	30% to 90% RH

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## **Customer contact information**

### **UK Customers only**

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd  
UK Customer Service & Support  
The Crescent  
Jays Close  
Basingstoke  
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: [micro.uksales@carefusion.com](mailto:micro.uksales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: [micro.ukservice@carefusion.com](mailto:micro.ukservice@carefusion.com)

Technical Support Enquiries:

Telephone: 01256 388551

Email: [support.rt.eu@carefusion.com](mailto:support.rt.eu@carefusion.com)

### **International customers only**

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH  
Customer Service & Support International  
Leibnizstrasse 7  
97204 Hoechberg  
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: [micro.internationalsales@carefusion.com](mailto:micro.internationalsales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: [support.admin.eu@carefusion.com](mailto:support.admin.eu@carefusion.com)

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: [support.rt.eu@carefusion.com](mailto:support.rt.eu@carefusion.com)

CareFusion UK 232 Ltd.,  
Quayside  
Chatham Maritime  
Kent ME4 4QY  
U.K.

**CE**  
0086

JDE no. 36-MAN1285  
Drg no. 048-14WW  
Issue 1.0  
February 2010  
© CareFusion 2010